

Traus

Ultrasonic Piezo



■ PIEZO TIP MANUAL

* Following some tips may be not included in the components.

Surgical operation	Using tip	Boost	Power	Pump	
Crestal approach	Sohn 028i Useful for sinus lift with water pressure		3	3	5
Sinus Lift	SB 400 Useful for sinus membrane separation inserting tip in the window of inside or outside		2	2	1
Window open	Sohn 016 Useful for opening window and operating osteotomy nearby mandible nerve.		3	3	3
Block bone	Sohn Saw, L/R Useful for Ramus Block bone collection And designed for Right and left type.	 L  R	3	3	3
Ridge expansion	Sohn Saw Useful for window open, block bone and ridge expansion		3	3	3
Implant site preparation / Bone harvest	SB 0320 Useful for bone collection, implant site preparation and alveolar bone osteoplasty		3	3	1
Extract	SB 0610 Useful for extraction of remained tooth and ridge expansion		3	3	3



The tip is not sterilized before packing. Please sterilize it before use.
Do not to let the Tip get to a patient's tongue, lips or cheeks.
It could cause a burn to patient.